







simcro ™
*Syringes,
Applicators
and
Accessories*





Animal Health
Delivery Systems


Y•TEX®

Simcro Syringes STF/STV

	DESCRIPTION	MAX. DOSES	CODE
 <p>1ml/2ml/5ml/10ml</p>	<p>Simcro STF Syringe</p> <ul style="list-style-type: none"> Extremely accurate repeating dose - variable dose adjustment. Durable construction, materials inert to virtually all actives and carriers. Award winning design. 	<p>1ml</p> <p>2ml</p> <p>5ml</p> <p>10ml</p>	<p>ST1</p> <p>ST2</p> <p>ST5</p> <p>ST10</p>
 <p>2ml/5ml</p>	<p>Simcro STF Syringe Bottle Mount</p> <ul style="list-style-type: none"> Extremely accurate repeating dose - variable dose adjustment. Durable, robust construction and materials. Suitable for 20 - 250ml glass bottles and plastic bottles. 	<p>2ml</p> <p>5ml</p>	<p>ST2B</p> <p>ST5B</p>
 <p>2ml/5ml/10ml 12ml/15ml</p>	<p>Simcro STV Syringe</p> <ul style="list-style-type: none"> Robust metal luer tip - no needle nut. Facilitates quick needle changes and eliminates breakages. Extremely accurate repeating dose. 12ml Draxxin syringe dose selected by cattle weight. High flow valve system which facilitates fast filling, especially on viscous products, is a standard feature on the 5ml, 10ml, 12ml and 15ml. 	<p>2ml</p> <p>5ml</p> <p>10ml</p> <p>12ml Draxxin</p> <p>12ml Hi Viscosity</p> <p>15ml</p>	<p>ST2V</p> <p>ST5V</p> <p>ST10V</p> <p>ST12V DX</p> <p>ST12V HV</p> <p>ST15V</p>
 <p>2ml/6ml</p>	<p>Simcro STV Syringe Bottle Mount</p> <ul style="list-style-type: none"> Extremely accurate repeating dose - simple dose adjustment. Durable, robust construction and materials. Suitable for 20 - 250ml glass bottles and plastic bottles. Available with bottle protecting shroud (SH). 	<p>2ml</p> <p>2ml both collars</p> <p>2ml shroud</p> <p>6ml</p> <p>6ml both collars</p> <p>6ml shroud</p>	<p>ST2V B</p> <p>ST2V B BC</p> <p>ST2V B SH</p> <p>ST6V B</p> <p>ST6V B BC</p> <p>ST6V B SH</p>
 <p>0.5ml 1.5ml</p>	<p>Simcro STV Syringe Selectable Dose Bottle Mount</p> <ul style="list-style-type: none"> Selectable dose - simple-to adjust-settings. Suitable for 20 - 250ml glass bottles and plastic bottles if purchased with both bottle collars (BC). Available in: 0.5ml, 5 stop, 0.1ml increments 1.5ml, 11 stop, 0.1ml increments 	<p>0.5ml</p> <p>0.5ml both collars</p> <p>0.5ml shroud</p> <p>1.5ml</p> <p>1.5ml both collars</p> <p>1.5ml shroud</p>	<p>ST0.5V B</p> <p>ST0.5V B BC</p> <p>ST0.5V B SH</p> <p>ST1.5V B</p> <p>ST1.5V B BC</p> <p>ST1.5V B SH</p>
	<p>Simcro STF Syringe & STV Replacement Barrels</p> <ul style="list-style-type: none"> Easy to replace. Clear dose adjustments marked. Replacement barrels are available for all types and sizes of syringes. 	<p>1ml STF</p> <p>2ml STF</p> <p>5ml STF</p> <p>10ml STF</p> <p>2ml STV</p> <p>5ml STV</p> <p>10ml STV</p> <p>2ml BM</p> <p>5ml STF BM</p> <p>2ml STV BM</p> <p>0.5ml STV BM</p> <p>1.5ml STV BM</p>	<p>ST1CYL</p> <p>ST2CYL</p> <p>ST5CYL</p> <p>ST10CYL</p> <p>ST2VMCYL</p> <p>ST5VMCYL</p> <p>ST10VMCYL</p> <p>ST2BCYL</p> <p>ST5BCYL</p> <p>ST2VBCYL</p> <p>ST0.5SVBCYL</p> <p>ST1.5SVBCYL</p>

Simcro Syringe Accessories


	DESCRIPTION	SIZE	CODE
	<p><i>Simcro Bottle Cradles</i></p> <ul style="list-style-type: none"> • To hold glass bottles or plastic containers securely around operator's neck or on livestock chute. • Complete with neck tie. • Leaves a hand free for easier operation. • Available in a range of sizes to fit 100ml, 250ml, and 500ml containers. 	100ml	NZBC100
		250ml	NZBC250
		500ml	NZBC500
 <p>20mm/27mm/33mm</p>	<p><i>Simcro Vented Draw-off</i></p> <ul style="list-style-type: none"> • To empty glass or plastic bottles where a vacuum restriction is likely to occur. • To fit bottles from 20ml to 500ml. 33mm comes with standard delivery tube. 20mm = 20-100 ml bottles 27mm = 250-500 ml US bottles 33mm = 250-500 ml European bottles 	20mm	D-O-S
		27mm	D-O-S
		33mm	D-O-XL
	<p><i>Upright Draw-off Spike</i></p> <ul style="list-style-type: none"> • Adapts to 100ml, 250ml & 500ml bottles. • Allows the bottle to remain upright in work tray or pocket. • Has built-in air inlet. 		EA102
 <p>60mm/100mm</p>	<p><i>Simcro Syringe Drench Nozzle</i></p> <ul style="list-style-type: none"> • Manufactured in nickle-plated brass tube with plastic nut. • For use with STF or special STV drench syringes. 60mm for goat and sheep. 100mm for horses. 	60mm (3 inches) 100mm (4 inches)	ST10NOZ S ST10NOZ



Simcro Flexi-Shot

- Enables user to inject livestock more remotely and enhances user safety as stock is kept at a distance.
- Can be used with most syringes.
- For intramuscular injection - especially swine.

CODE ST FLEXISHOT




Simcro Handi-Ject

- Allows for better control and visibility while injecting.
- Allows subcutaneous and intramuscular injection.
- Can be used with most syringes.
- Enables user to inject stock more remotely avoiding potential injury.

CODE ST HANDIJECT

Simcro Horse Drenchers



Simcro STV Horse Drencher

- Specialized equine drencher.
- 100mm drench nozzle.
- Durable plated brass construction.
- Can be fitted with different size bottles/containers.

10ml

MAX. DOSE **CODE**

10ml ST10V D



Simcro STF Horse Drencher

- Specialized equine drencher.
- 100mm drench nozzle.
- Durable plated brass construction.
- Available with 125ml bottles.

10ml



MAX. DOSE **CODE**

10ml ST10 EB


Simcro Applicators Pour-on/Drenchers

	DESCRIPTION	MAX. DOSES	CODE	
15ml/30ml	Simcro Applicator Pour-on Gun <ul style="list-style-type: none"> Extremely reliable and durable. A highly respected and sought-after applicator. Suitable for sheep and goats. Available with 6.4mm and 9.5mm hose barb. 	15ml	ST15 P	
		30ml	ST30 P	
60ml	Simcro Applicator Pour-on Gun <ul style="list-style-type: none"> New rose tip nozzle reduces wastage and post delivery spurt by 90%. Large inlet Hi Flow system allows fast barrel refilling. Handles viscous materials in cold conditions. Available with a 9.5mm hose barb. 	60ml	ST60 P	
15ml/30ml	Simcro Applicator Drencher Gun <ul style="list-style-type: none"> Accurate and reliable delivery system. Strong durable construction. Both nozzles are manufactured in durable plated brass (100mm-sheep, 150mm-cattle). Available with 6.4mm and 9.5mm hose barb. The 15ml Simcro Drencher is an industry leader. 	NOZZLE		
		15ml	100mm	ST15 DS
		15ml	150mm	ST15 D
		30ml	100mm	ST30 DS
		30ml	150mm	ST30 D

Simcro Applicator Accessories

	Simcro Pour-on Nozzle <ul style="list-style-type: none"> Fits 15ml and 30ml units. Single stream applicator. 		Simcro Fan Spray Pour-on Nozzle <ul style="list-style-type: none"> Fits 15ml and 30ml units. Spreads out pour-on over animal's back. Durable plated brass construction.
	CODE ST25NOZ P		CODE ST25NOZ F
100mm	Simcro Calf/Sheep Drench Nozzle <ul style="list-style-type: none"> Fits 15ml and 30ml units. Use for calf, sheep and goats. Durable plated brass construction. 100ml length. 	150mm	Simcro Cattle Drench Nozzle <ul style="list-style-type: none"> Fits 15ml and 30ml units. Use for cattle, 150mm length. Durable plated brass construction.
	CODE ST25NOZ S		CODE ST25NOZ C

Simcro Backpacks

	DESCRIPTION	SIZE	CODE
	<ul style="list-style-type: none"> Free standing for safe, easy filling and display purposes. 2.5L & 5L backpacks provide a unique cavity for spigot cap and strap. Packs allow for overfilling to cater for promotions where additional product is supplied. Comes complete with straps, transit cap and spigot cap for inverted use on back or chute. 	1 Liter	STBP10
		2.5 Liter	STBP25
		5 Liter	STBP50

1L/2.5L/5L*

*L=Liter (3.8 liters = 1 Gallon)

Simcro Backpack Accessories



DESCRIPTION

Simcro Upright Draw-off

- Suitable for any anthelmintic, pour-on or spray product.
- Customized to fit any cap or container.
- Vented to allow easy evacuation of the product from backpack, to avoid collapsing of backpack and also to avoid vacuum build-up which makes application inaccurate, slow and difficult.

SPIGOT SIZE	BACKPACK SIZE	CODE
6.4mm	2.5 Liter	BP CAP U 25
6.4mm	5.0 liter	BP CAP U 50
	GALLON CONTAINER	
6.4mm	Up to 1 gallon	NZ CAP GAL

Simcro 38mm Inverted Spigot Caps



- Designed for use with inverted backpacks to allow liquid to be drawn from the bottom of the container.
- Available in two stage (6.4-9.5mm) to fit onto tubing.
- Vented to allow easy evacuation of the product from backpack, to avoid collapsing of backpack and also to avoid vacuum build-up which makes application inaccurate, slow and difficult.

SPIGOTS	CODE
---------	------

6.4 mm/9.5mm	BPCAP STEP
9.5mm	BPCAP WIDE


Simcro Backpack Straps



- To be threaded on backpack, enabling unit to be worn on user's back.
- Also can be used to hang near the chute head for holding a container (1L, 2.5L, 5L).
- Strong and durable. Webbed tape 2500mm long x 25mm wide incorporating a plastic buckle.

CODE	BPSTRAP
------	---------

Simcro 38mm Cap Transit



- To be used with Simcro Backpacks.

SIZE	CODE
------	------

38mm	BPCAP TRAN
------	------------

Simcro Tube Adaptor





























- To fit onto tubing to allow the connection between various internal diameters of tube.
- Available in size: 4.6mm - 6.4mm

SIZE	CODE
------	------

4.6mm - 6.4mm	ST JOIN 1
---------------	-----------

Simcro Combinations

STF								
STV								
Accessories								
Key: 1. Small draw-off fits 50 & 100ml bottles (inverted) 2. Large draw-off fits 250 & 500ml bottles (inverted) 3. Extra Large draw-off fits 250 & 500ml European bottles 4. Spike draw-off for use with all sized Upright bottles								
STF	1ml	ST1	ST1S	ST1L	ST1SL		ST1U	
	2ml	ST2	ST2 S	ST2 L	ST2 SL	ST2 XL	ST2 U	ST2 B
	5ml	ST5	ST5 S	ST5 L	ST5 SL	ST5 XL	ST5 U	ST5 B
	10ml	ST10	ST10 S	ST10 L	ST10 SL	ST10 XL	ST10 U	
STV	0.5ml							ST0.5SV B
	1.5ml							ST1.5SV B
	2ml	ST2V	ST2V S	ST2V L	ST2V SL	ST2V XL	ST2V U	ST2V B
	5ml	ST5V	ST5V S	ST5V L	ST5V SL	ST5V XL	ST5V U	
	6ml							ST6V B
	10ml	ST10V	ST10V S	ST10V L	ST10V SL	ST10V XL	ST10V U	
	12ml				ST12SVHV SL			
	15ml	ST15V	ST15V S	ST15V L	ST15V SL	ST15SV XL	ST15V U	

Manufactured by:
 Simcro Tech Limited
 6-8 Latham Court, Box 5152,
 Hamilton, New Zealand
 Phone: +64 7 846 7723 Fax: +64 7 846 7738
 Web: www.simcro.com
 E-Mail: info@simcro.com



Distributed by:
 Y-TEX Corporation
 1825 Big Horn Avenue,
 Cody, WY 82414
 Phone: 888-600-YTEX
 Web: www.ytex.com
 E-Mail: ytexinfo@ytex.com